

Plant and Animal Cells Lab Activity



Life Science
**Plant and
Animal Cells**
Lab Activity



Plant Cells vs Animal Cells Lab

Activity Learning Objectives:

At the end of this laboratory activity, students are expected to:

- Compare the structure of plant cells (onion) and animal cells (cheek cell).
- Cite the similarities and differences between plants and animal cells.

Plant Cells vs Animal Cells Lab


Activity Materials:

Materials:

- Microscope
- Coverslips
- Slides
- Lens paper
- Onion
- Cheek cells
- Elodea – an aquatic plant found in aquariums
- Water
- Dropper
- Iodine
- Methylene blue
- Toothpicks
- Forceps
- Colored pencils


2. For the onion:

- Remove a scale from the onion.
- Peel a thin layer of onion from its inner surface.
- Make a wet-mount slide of the onion tissue, then add 2-3 drops of water.
- Put on a coverslip, then stain with iodine.
- Observe, draw what you can see, and label what you can recognize.





KINEMASTER







KINEMAS




4. For the cheek:

- Gently scrape the inside of your mouth using a flat clean toothpick. Don't worry if you can't see anything. You can actually pick up any cheek cell when you scrape.
- Smear the toothpick on the slide to spread the saliva and cells. You may add a tiny drop of water if the sample is too dry.
- Put on a coverslip, and stain with methylene blue.
- Observe, draw what you can see, and label what you can recognize.







Animal cells under microscope
at 10x (low power)



4:26 ✓✓



Animal cells under microscope at
100x(high power)

4:33 ✓✓

Name _____ Date _____



Find the function of all the labeled parts of these cells and to your notebook.

